## acher.co.ke

Download this and other FREE materials from https://teacher.co.ke/notes



## Each candidate should have the following

- 20Ml of solution Q
- 10% sodium hydroxide solution
- 1% copper II sulpahte solution
- 0.1% DCPIP
- Filter paper
- 4 test tubes in a rack
- Source of heat
- 1 dropper

## Preparation

Solution Q is a solution mixture of egg albumen, Olive oil and ascorbic acid in the ratio of 1:1:1 respectively.

This paper consists of 1 printed page