Name	Adm No
Class	
233/3	Date
CHEMISTRY	
Paper 3	1 CO COLLENA C
JUNE 2022	MARKING SCHEME

Time: 21/4 Hours

# SUKELLEMO FORM 4 CHEMISTRY PRACTICAL

# **INSTRUCTIONS**

- \* Write your name and the indicated details in the spaces provided on this page.
- \* Answer ALL questions in the spaces provided in the question paper.
- \* You are **NOT** allowed to start working with the apparatus for the *first 15 minutes* of the 2<sup>1</sup>/<sub>4</sub> hours allowed for this paper. This time is to enable you to read the question paper and make sure you have all the chemicals and apparatus you may need.
- \* Mathematical tables and silent electronic calculators may be used.
- \* All working must be clearly shown where necessary.

# For Examiner's Use Only

Question	Maximum score	Candidate's score
1	21	
2	10	
3	9	
Total score	40	

This paper has 8 printed pages. Candidates should check the question paper to ascertain that <u>ALL</u> pages are printed as indicated and that no questions are missing.

# **QUESTION 1**

# You are provided with:

- 6.0g of an alkanoic acid labelled solid P in a boiling tube.
- 2M sodium hydroxide solution labelled solution Q.

# You are required to:

- (a) Determine the solubility of solid P at different temperatures.
- (b) Determine the number of moles of water of crystallization in solid P.
- (c) Find the molar mass of the alkanoic acid.

#### Procedure I

- i) Using a burette add 10cm<sup>3</sup> of distilled water to the solid P in the boiling tube. Heat the mixture while stirring with a thermometer to about 70°C. When the entire solid has dissolved allow the solution to cool while stirring with the thermometer. Note the temperature at which crystals of solid P first appear. Record this temperature in the Table 1 below.
- ii) Using the burette add 2cm<sup>3</sup> of distilled water to the contents of the boiling tube. Warm the mixture, while stirring with the thermometer until all the solid dissolves. Allow the mixture to cool while stirring. Note the temperature at which the crystals of solid P first appear.
- iii) Repeat the procedure (ii) two more times and record the temperatures in Table I. Retain the contents of the boiling tube for use in procedure II.

# TABLE I

Volume of water in the boiling tube (cm <sup>3</sup> )	Temperature at which crystals of solid P first appear (°C)	Solubility of Solid P (g/100gH <sub>2</sub> O)
10	55	6×100 = 60
12	50	$\frac{6\times 100}{12} = 50$
14	45	6×100 = 42.9
16	40	6×100 = 37.5

a) i) Complete Table 1 by calculating the solubility of solid P at different temperatures (6 mks)

- Calculation


- School value \$ 2°C - Im K (whose school value is 1st reading) - Ink

- Trend (decrease in temperature)






ii) On the grid provided, plot a graph of the solubility of solid P (vertical axis) against temperature. (3 mks)



iii) Using the graph determine the temperature at which 55g of solid P would dissolve in 100cm<sup>3</sup> of water. (1 mk)



Page 3 of 8

# Procedure II

i) Transfer the contents of the boiling tube in Procedure I into a 250 ml volumetric flask. Rinse both the boiling tube and the thermometer with distilled water and add to the volumetric flask. Add more distilled water to make up to the mark. Transfer the solution to a 250 ml beaker, Label this solution R. Wash the volumetric flask then rinse it with distilled water ready for use in Step(ii).

ii) Using a measuring cylinder place 25.0 cm<sup>3</sup> of solution Q into a 250 ml volumetric flask. Add about 200cm<sup>3</sup> of distilled water. Shake well. Add more distilled water to the make up to the mark. Label this solution T.

iii) Fill a burette with solution R. Using a pipette and pipette filler, place 25 cm<sup>3</sup> of solution T into a 250 ml conical flask. Add 2-3 drops of phenolphthalein indicator and titrate with solution R. Record your results in table II. Repeat the titration two more times and complete the table.

# TABLE II

Titre	1	2	3
Final burette reading (cm <sup>3</sup> )	13.2	27.2	13.2
Initial burette reading (cm <sup>3</sup> )	0:0	14,0	0.0
Volume of solution R used (cm <sup>3</sup> )	13.2	13.2	13.2

(4 marks)

DETERMINE:

Values Averaged pust by profitin 20.2 of

a) The average volume of solution R used (1 mk)

b) The concentration of solution T in moles per litre. (1 mk)

1000 cm² => 2 moles 0.05 moles => 250 cm² => 1000 cm² => 1000 cm²

1000 x 0,02 = 0.2 H

Page 4 of 8

6

c) The concentration of alkanoic aci-	d solution R in r	noles per litre (1 mole	of acid reacts with 2
moles of the base)			
Holory bare in scale > (D)	Acid	(2 mks)	$\sim$
	0.0025	values => Av	Corrections
1000 = 0.005 rules	A VALUE OF THE PROPERTY OF THE	Voline	water by
M.O. (5th) (a)		2	Penalite
Mikraho of Addibase	e'e	10000	Penalite  -1/2 Mt pur worn a  or mussing un
	- 1000	cm3 x0.002	) or mussing un
0.00x 1 0.00sng	Au	volume V	
2 = 0.000 rushis air			
d) The relative formula mass of the	alkanoic acid, so	olid P (2 mks)	
		)	
Ans in (1) above x	350	L correct au	n
1422			ectous.
Aus in Coabre x 250 0	<u> </u>	de curit	CITOUS,
e) The formula of P has the form M	LxH <sub>2</sub> O Determ	ine the value of x in the	e formula. (The
relative formula mass of M is 90; C		A CONTRACT OF THE PROPERTY OF	
RFM of Pin	H. 2 H20 :	- Am in (d)	above.
	x = Ans	in(d) above	-90ND
		1.8	The same of the sa
			10
	= Co	red aus. L	
		Dr. 12 Pill 3 12 1 V.	TA VIEW NEW AND
Conclize -	- LMK 11	the value of	ru a
Penalize - Jactus	2	TO CAR TO A DECEMBER OF THE PARTY OF THE PAR	a coleibas
Tracerus			
(N. V.6) - 18		Kirk Company	The state of the s
		Kurinen tops	
	TYPE K	SPINE DESCENDE	revolute at the

# **QUESTION 2**

You are provided with solid Y. Use it to carry out the tests indicated below and record your observations and inferences in the spaces provided.

a) Put half a spatula of solid Y in a boiling tube. Add about 5 cm<sup>3</sup> of distilled water and shake. Divide the resultant mixture into four portions of 1 cm<sup>3</sup> each.

Observations	Inferences
Solid dissolves to form	- Soluble Sall.
a Colonvless Solvikon (1 mk)	Colored contakes (1 mk)
(1 mk)	Cu2+ Fe2+ Fe3+ (1 mk)

i) To the first portion add aqueous sodium hydroxide dropwise till in excess.

Observations	Inferences
white precipitate dissolves in expens) to form a colourless solution (1 mk)	Zn2+ A13+ pt2+  () present (1 mk)

ii) to the second portion add aqueous ammonia dropwise till in excess

Observations	Inferences
white precipitate dissolves in excess to popular colourless solution (1 mk)	Zn2+ present (1) (1 mk)

iii) To the third portion add lead (II) nitrate solution and warm.

Observations	Inferences
not dissolve in warming	> SO42-, SO32-, CO32-100 sent
The getter processing (I'm)	
The house backling of the	(1 absent (+mk)
	tied to correct observation

Page 6 of 8



iv) To the fourth portion add dilute nitric (V) acid followed by Barium nitrate solution.

Observations	Inferences
- No enewercencer to	SO4 present
(1 mk)	(1 mk)

# **QUESTION 3**

You are provided with solid M and solid F.

Carry out the tests below and write the observations and inference in the spaces below:

a) Using a metallic spatula, ignite about one half of solid M in a Bunsen burner flame.

Observations	Inferences
Solid melts and burns with a yellow sooty Home.	C=C - C=C
(1 mk)	(1 mk)

b) Place the other half of solid M in a boiling tube. Add 15cm<sup>3</sup> of distilled water and shake well. Label this solution M. Use this solution for the following tests.

i) Place 2cm<sup>3</sup> of solution M in test tube and determine pH.

Observations	Inferences
pH=3	Strongly acidic M
ii) To about 2 cm <sup>3</sup> of the solution M obtained in	(½ mk)

manganate (VII)

Observations	Inferences
Purple KHOG (ap) thoughs	C=C! . C=C
	prosent ROH
(1 mk)	(1 mk)
	I A would be private to be



iii) To about 2 cm<sup>3</sup> of the solution M obtained in (b) above add 2 drops of bromine water. and

Observations	Inferences
Yellow browne water changes colourless (1/2 mk)	Present (1/2 mk)

iv) To the remaining solution M in the boiling tube, add half of solid N-F

Observations		Inferences
Effervence		R-CZO PCCEPT RCOOH  Present (10)  - COOH
	(½ mk)	(½ mk)